

SERPINA10 Rabbit pAb

货号: B16470

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100
理论分子量	50kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse liver,Rat liver
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 165-444 of human SERPINA10 (NP_001094077.1).
序列	FAFIHKDFDVKETFFNLKRYFDTECVPMNFRNASQAKRLMNHYINKETRGKIPKLFDEINPETKLILVDYILFKGKWLT PFD PVFTEVDTFHLDKYKTIKVPMMYGAGKFASTFDKNFRCHVLKLPYQGNATMLVVLMEKMGDHLALEDYLTDDL VETWLRN MKTRNMEVFFPKFKLDQKYEMHELLRQMGIRRI FSPFADLSELSATGRNLQVSRVLQRTVIEVDERGTEAVAGILSEITAYS MPPVIKVD RPFHFMIYEETSGMLLFLGRVVNPTLL

靶点信息

研究背景	The protein encoded by this gene belongs to the serpin family. It is predominantly expressed in the liver and secreted in plasma. It inhibits the activity of coagulation factors Xa and XIa in the presence of protein Z, calcium and phospholipid. Mutations in this gene are associated with venous thrombosis. Alternatively spliced transcript variants have been found for this gene.
------	---

基因ID	51156
基因名	SERPINA10
Swiss	Q9UK55
别名	SERPINA10;PZI;ZPI

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn