

LDLRAP1 Rabbit pAb

货号: B16476

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:100
理论分子量	33kDa
实测分子量	36KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,Mouse liver
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 59-308 of human LDLRAP1 (NP_056442.2).
序列	EQPKGEELSAAAIAKRIVATAKASGKKLQKVTLKVSPRGILTDNLTNQLIENVSIYRISYCTADKMHDKVFAYIAQSQHNQSLECHAFLECTKRKMAQAVTLTVAQAFKVAFEFWQVSKEEKEKRDKASQEGGDVLGARQDCTPSLKSLVATGNLLDLEETAKA PLSTVSANTTNMDEVPRPQALSGSSVWELDDGLDEAFSRLAQSRITNPQVLDTGTLTAQDMHYAQCLSPVDWDKPDSS GTEQDDLFSF

靶点信息

研究背景	The protein encoded by this gene is a cytosolic protein which contains a phosphotyrosine binding (PTD) domain. The PTD domain has been found to interact with the cytoplasmic tail of the LDL receptor. Mutations in this gene lead to LDL receptor malfunction and cause the disorder autosomal recessive hypercholesterolemia.
------	--

基因ID	26119
基因名	LDLRAP1
Swiss	Q5SW96
别名	LDLRAP1;ARH;ARH1;ARH2;FHCB1;FHCB2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn