

L3MBTL1 Rabbit pAb

货号: B16479

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:100
理论分子量	43kDa/82kDa/83kDa/85kDa/92kDa
实测分子量	43kDa/82kDa/83kDa/85kDa/92kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 533-772 of human L3MBTL1 (NP_056293.4).
序列	PLSYRSLPHTRTSKYSFHHRKCPTPGCDGSGHVTGKFTAHHCLSGCPLAERNQSRLKAELSDSEASARKKNLSGFSPRKKPRHHGRIGRPPKYRKIPQEDFQTLTPDVVHQSLFMSALSAHPDRSLVCWEQHCKLLPGVAGISASTVAKWTIDEVFGFVQTLTGCEDQARLFKDEARIVRVTHVSGKTLVWTVVAQLGDLVCSDDLQEGKGILETGVHSLCPLTHLLAKLSFASDSQY

靶点信息

研究背景	This gene represents a polycomb group gene. The encoded protein functions to regulate gene activity, likely via chromatin modification. The encoded protein may also be necessary for mitosis. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	26013

基因名	L3MBTL1
Swiss	Q9Y468
别名	L3MBTL1;H-L(3)MBT;L3MBTL;ZC2HC3;dj138B7.3

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn