

RRAS2 Rabbit pAb

货号: B16484

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	15kDa/19kDa/23kDa/24kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,SKOV3,NIH/3T3,Mouse kidney,Mouse heart,Mouse lung,Rat kidney
细胞定位	Cell membrane,Cytoplasmic side,Lipid-anchor
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-204A of human RRAS 2 (NP_036382.2).
序列	MAAAGWRDGGSGQEKYRLVVVGGGGVGKSALTIQFIQSYFVTDYDPTIEDSYTKQCVIDDRAARLDILDTAGQEFGAMR EQYMRTGEGFLLVFSVTDGRGSFEEIYKFQRQILRVKDRDEFPMILIGNKADLDHQRQVTQEEGQQLARQLKVTYMEASAKI RMNVDQAFHELVRVIRKFQEQECPPSPEPTRKEKDKKGCHCVIF

靶点信息

研究背景	This gene encodes a member of the R-Ras subfamily of Ras-like small GTPases. The encoded protein associates with the plasma membrane and may function as a signal transducer. This protein may play an important role in activating signal transduction pathways that control cell proliferation. Mutations in this gene are associated with the growth of certain tumors. Pseudogenes of this gene are found on chromosomes 1 and 2. Alternate splicing results in multiple transcript variants.
------	---

基因ID	22800
基因名	RRAS2
Swiss	P62070
别名	RRAS2;TC21

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn