

NAA10 Rabbit pAb

货号: B16514

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:100
理论分子量	24kDa/26kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,HeLa,HL-60,HepG2,HT-29,Mouse thymus,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 161-235 of human NAA 10 (NP_003482.1).
序列	HLELKEKGRHVVLGAIENKVESKGNPPSSGEACREEKGGLAAEDSGGDSKDLSEVSETTESTDVKDSSEASDSAS

靶点信息

研究背景	N-alpha-acetylation is among the most common post-translational protein modifications in eukaryotic cell s. This process involves the transfer of an acetyl group from acetyl-coenzyme A to the alpha-amino group on a nascent polypeptide and is essential for normal cell function. This gene encodes an N-terminal acetyltransferase that functions as the catalytic subunit of the major amino-terminal acetyltransferase A complex. Mutations in this gene are the cause of Ogden syndrome. Alternate splicing results in multiple transcript variants.
------	--

基因ID	8260
基因名	NAA10
Swiss	P41227
别名	NAA10;ARD1;ARD1A;ARD1P;DXS707;MCOPS1;NATD;OGDNS;TE2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn