

UBE2V2 Rabbit pAb

货号: B16519

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	16kDa
实测分子量	16kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,SKOV3,HeLa,HL-60,Mouse liver,Mouse skin,Rat liver,Rat brain
细胞定位	cytoplasm,extracellular exosome,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-145 of human UBE2V 2 (NP_003341.1).
序列	MAVSTGVKVPNFRLLLEELEGQKGVGDGTVSWGLEDDEDMTLTRWTGMIIGPPRTNYENRIYSLKVECGPKYPEAPPSV RFVTKINMNGINNSSGMVDARSIPVLAKWQNSYSIKVVLQELRRLMMSKENMKLPQPPEGQTYNN

靶点信息

研究背景	Ubiquitin-conjugating enzyme E2 variant proteins constitute a distinct subfamily within the E2 protein family. They have sequence similarity to other ubiquitin-conjugating enzymes but lack the conserved cysteine residue that is critical for the catalytic activity of E2s. The protein encoded by this gene also shares homology with ubiquitin-conjugating enzyme E2 variant 1 and yeast MMS2 gene product. It may be involved in the differentiation of monocytes and enterocytes.
------	---

基因ID	7336
基因名	UBE2V2
Swiss	Q15819
别名	UBE2V2;DDVIT1;DDVit-1;EDAF-1;EDPF-1;EDPF1;MMS2;UEV-2;UEV2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn