

BMP3 Rabbit pAb

货号: B16571

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:100 - 1:200
理论分子量	53kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,Mouse ovary,Mouse lung
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 243-472 of human BMP 3 (NP_001192.2).
序列	YANDAAISEPESVSSSLQGHRNFPTGTVPKWDSHIRAALSIERRKKRSTGVLLPLQNNELPGAHEYQYKKDEVWEERKPYKT LQAQAPEKSKNKKKQKRGPHRKSQLQFDEQTLKKARRKQWIEPRNCARRYLKVDFADIGWSEWIISPKSFDAYYCSGAC QFPMPKSLKPSNHATIQSIVRAVGVPPIPEPCCVPEKMSSLSILFFDENKNVVLKVYPNMTVESCACR

靶点信息

研究背景	BMP3 belongs to the transforming growth factor-beta (TGFB) superfamily. Bone morphogenic protein, als o known as osteogenin, induces bone formation.
基因ID	651

基因名	BMP3
Swiss	P12645
别名	BMP3;BMP-3A

产品验证



Western blot analysis of BMP3 expressed in SW480,Mouse ovary,Mouse lung using BMP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn