

UBE2T Rabbit pAb

货号: B16581

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	22kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse thymus,Rat spleen,Rat heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-197 of human UBE2T (NP_054895.1).
序列	MQRASRLKRELHMLATEPPPGITCWQDKDQMDDLRAQILGGANTPYEKGVFKLEVIIPERYPFEPPIRFLTPIYHPNIDSAG RICLDVLKLPPKGAWRPSLNIATVLTSLQLLMSEPNPDDPLMADISSEFKYNKPAFLKNARQWTEKHARQKQKADEEEMLD NLPEAGDSRVHNSTQKRKASQLVGIEKKFHPDV

靶点信息

研究背景	The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. D effects in this gene have been associated with Fanconi anemia of complementation group T. Two transcrip t variants encoding different isoforms have been found for this gene.
基因ID	29089

基因名	UBE2T
Swiss	Q9NPD8
别名	UBE2T;FANCT;HSPC150;PIG50

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn