

CCS Rabbit pAb

货号: B16587

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	29kDa
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human CCS (N P_005116.1).
序列	MASDSGNQGTLCLEFAVQMTCQSCVDAVRKSLQGvagvqdvevhLEDQMVLVHTTLPSQEVQALLEGTRQAVLKG MSGGQLQNLGA AVAILGGPGTVQGVVRFQLTPERCLIEGTIDGLEPGLHGLHVHQQYGDLTNNCNSCGNHFNPDGASHG GPQSDRHRGDLGNVRADADGRAIFRMEDEQLKVWDVIGRSLIIDEGEDDLGRGGHPLSKITGNSGERLACGIIARSAGL FQNPQKICSCDGLTIWEERGRPIAGKGRKESAQPPAHL

靶点信息

研究背景	Copper chaperone for superoxide dismutase specifically delivers Cu to copper/zinc superoxide dismutase and may activate copper/zinc superoxide dismutase through direct insertion of the Cu cofactor.
基因ID	9973

基因名	CCS
Swiss	O14618
别名	CCS

产品验证



Western blot analysis of CCS expressed in 293T using CCS Rabbit pAb at 1:1000. Secondary antibody: HR P Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn