

USP1 Rabbit pAb

货号: B16593

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	88kDa
实测分子量	115kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HepG2,MCF7,Mouse testis,Mouse spleen,Mouse brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human USP1 (NP_003359.3).
序列	MPGVIPSESNGLSRGSPSKNRLSLKFFQKKETKRALDFTDSQENEEKASEYRASEIDQVVPAAQSSPINCEKRENLLPFVG LNNLGNTCYLNSILQVLYFCPGFKSGVKHLFNIIISRKKEALKDEANQKDKGNCKEDSLASYELICSLQSLIISVEQLQASFLN PEKYTDELATQPRRLNLTRELNPMYEGYLQHDAQEVLQCILGNIQETCQLLKKEEVKNVAELPTKVVEIIPHPKEEMNGINSIE MDSMRHSEDFKEKLPKGNGKRKSDTEFGNM

靶点信息

研究背景	This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases that is a deubi quitinating enzyme (DUB) with His and Cys domains. This protein is located in the cytoplasm and cleaves the ubiquitin moiety from ubiquitin-fused precursors and ubiquitinylated proteins. The protein specifically deubiquitinates a protein in the Fanconi anemia (FA) DNA repair pathway. Alternate transcriptional splice variants have been characterized.
------	---

基因ID	7398
基因名	USP1
Swiss	O94782
别名	USP1;UBP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn