

SPINT2 Rabbit pAb

货号: B16609

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	21kDa/28kDa
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 28-197 of human SPINT 2 (NP_066925.1).
序列	ADRERSIHDFCLVSKVVGRCRASMPRWYNVTDGSCQLFVYGGCDGNSNNYLTKEECLKKCATVTENATGDLATSRNA ADSSVPSAPRRQDSEDHSSDMFNYYEYCTANAVTGPCRASFPRWYFDVERNSCNNFIYGGCRGNKNSYRSEEACMLRC FRQQENPPLPLGSK

靶点信息

研究背景	This gene encodes a transmembrane protein with two extracellular Kunitz domains that inhibits a variety of serine proteases. The protein inhibits HGF activator which prevents the formation of active hepatocyte growth factor. This gene is a putative tumor suppressor, and mutations in this gene result in congenital so dium diarrhea. Multiple transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	10653
基因名	SPINT2
Swiss	O43291
别名	SPINT2;DIAR3;HAI-2;HAI2;Kop;PB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn