

MYLK3 Rabbit pAb

货号: B16638

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	52kDa/88kDa
实测分子量	52kDa/88kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human MYLK 3 (NP_872299.2).
序列	MQQDAAQHGARLEALFRMVAAVDRAIALVGATFQKSKVADFLMQGRVPWRRGSPGDSPEENKERVEEEGGKPKHVLSTSGVQSDAREPGEESQKADVLEGT

靶点信息

研究背景	Phosphorylation of cardiac myosin heavy chains (see MYH7B, MIM 609928) and light chains (see MYL2, MIM 160781) by a kinase, such as MYLK3, potentiates the force and rate of cross-bridge recruitment in cardiac myocytes (Chan et al., 2008 [PubMed 18202317]).[supplied by OMIM, Jul 2008]
基因ID	91807

基因名	MYLK3
Swiss	Q32MK0
别名	MYLK3;MLCK;MLCK2;caMLCK

产品验证



Western blot analysis of MYLK3 expressed in CACO2,mouse brain using MYLK3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn