

MMRN1 Rabbit pAb

货号: B16642

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	58kDa/138kDa
实测分子量	180kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,U-251MG,Mouse lung,Mouse spleen
细胞定位	Secreted
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 949-1228 of human MMRN1 (NP_031377.2).
序列	IKHSLPDIQLLQKGLTEFVEPIIQIKTQAALSNLTCCIDRSLPGSLANVVKSQKQVKSLPKKINALKKPTVNLTTVLIGRTQRNT DNIIYPEEYSSCSRHPCQNGGTCINGRTSFTCACRHPFTGDNCTIKLVEENALAPDFSKGSYRYAPMVAFFASHTYGMTIPG PILFNNLDVNYGASYTPRTGKFRIPYLGYYVFKYTIESFSAHISGFLVVDGIDKLAFESENINSEIHCDRVLTGDALLELNYGQ EVLRLAKGTIPAKFPPVTTFSGYLLYRT

靶点信息

研究背景	Multimerin is a massive, soluble protein found in platelets and in the endothelium of blood vessels. It is comprised of subunits linked by interchain disulfide bonds to form large, variably sized homomultimers. Multimerin is a factor V/Va-binding protein and may function as a carrier protein for platelet factor V. It may also have functions as an extracellular matrix or adhesive protein. Recently, patients with an unusual autosomal-dominant bleeding disorder (factor V Quebec) were found to have a deficiency of platelet multimerin.
基因ID	22915
基因名	MMRN1
Swiss	Q13201
别名	MMRN1;ECM;EMILIN4;GPIa*;MMRN

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn