

MFI2 Rabbit pAb

货号: B16645

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	32kDa/80kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,SW480,A375,Mouse liver,Mouse brain,Mouse thymus
细胞定位	Cell membrane,GPI-anchor,Lipid-anchor
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of human MFI2 (NP_201573.1).
序列	GMEVRWCATSDPEQHKCGNMSEAFREAGIQPSLLCVRGTSADHCVQLIAAQEADAITLDGGAIYEAGKEHGLKPVVGEVY DQEVGTSYYAVAVVRRSSHVTIDTLKGVKSCHTGINRTVGWNVPGYLVESGRLSVMGCDVLKAVSDYFGGSCVPGAG ETSYSESLCRLCRGDSSGEGVCD

靶点信息

研究背景	The protein encoded by this gene is a cell-surface glycoprotein found on melanoma cells. The protein shares sequence similarity and iron-binding properties with members of the transferrin superfamily. The importance of the iron binding function has not yet been identified. This gene resides in the same region of chromosome 3 as members of the transferrin superfamily. Alternative splicing results in two transcript variants.
------	--

基因ID	4241
基因名	MELTF
Swiss	P08582
别名	MELTF;CD228;MAP97;MFI2;MTF1;MTf

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn