

CMPK1 Rabbit pAb

货号: **B16681**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:500 - 1:1000
理论分子量	16kDa/22kDa
实测分子量	22kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HL-60,HT-29,HeLa,Mouse brain,Mouse kidney
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 69-228 of human CMPK 1 (NP_057392.1).
序列	LLRDERKNPDSQYGELIEKYIKEGKIVPEITISLLKREMDQTMANAQKNKFLIDGFPRNQDNLQGWNKTM MDGKADV SFV LFFDCNNEICIERCLER GKSSGR SDDNRESLEKRIQTYLQSTKPIIDLYEEMGKVKKIDASKSVDEVVQIFDKEG

靶点信息

研究背景	This gene encodes one of the enzymes required for cellular nucleic acid biosynthesis. This enzyme catalyzes the transfer of a phosphate group from ATP to CMP, UMP, or dCMP, to form the corresponding diphosphate nucleotide. Alternate splicing results in both coding and non-coding transcript variants.
基因ID	51727

基因名	CMPK1
Swiss	P30085
别名	CMPK1;CK;CMK;CMPK;UMK;UMP-CMPK;UMPK

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn