

CHIA Rabbit pAb

货号: B16686

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	33kDa/40kDa/52kDa
实测分子量	50-55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HepG2,Mouse stomach,Rat lung,Rat stomach
细胞定位	Cytoplasm,Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-368 of human CHIA (N P_001244932.1).
序列	MVSTPENRQTFITSVIKFLRQYEFDGLDFDWEYPGSRGSPQDKHLFTVLVQEMREAFEQEAKQINKPRLMVTAAVAAGIS NIQSGYEIPQLSQYLDYIHVMYDYLHGSWEGYTGENSPLYKYPTDTGSNAYLNVDYVMNYWKDNGAPAEKLIVGFPTYGH NFILSNPSNTGIGAPTSAGAPAGPYAKESGIWAYYEICTFLKNGATQGWDAPQEVYAYQGNVWVGYDNIKSFDIKAQWL KHNFKGAMVWAILDDFTGTFCNQGKFPLISTLKKALGLQSASCTAPAQPIEPITAAPSGSGNGSGSSSSGGSSGGSGF CAVRANGLYPVANNRNAFWHCVNGVTYQQNCQAGLVFDTSCDCCNWA

靶点信息

研究背景	The protein encoded by this gene degrades chitin, which is found in the cell wall of most fungi as well as in arthropods and some nematodes. The encoded protein can also stimulate interleukin 13 expression, and variations in this gene can lead to asthma susceptibility. Several transcript variants encoding a few different isoforms have been found for this gene.
------	--

基因ID	27159
基因名	CHIA
Swiss	Q9BZP6
别名	CHIA;AMCASE;CHIT2;TSA1902;chitinase;acidic

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn