

AGPAT2 Rabbit pAb

货号: B16698

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	27kDa/30kDa
实测分子量	31kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29
细胞定位	Endoplasmic reticulum membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 199-278 of human AGPAT2 (NP_006403.2).
序列	IVPVVYSSFSSFYNTKKKFFTSGLTVTVQVLEAIPTSGTLTAADVPALVDTCHRAMRTTFLHISKTPQENGATAGSGVQPAQ

靶点信息

研究背景	This gene encodes a member of the 1-acylglycerol-3-phosphate O-acyltransferase family. The protein is located within the endoplasmic reticulum membrane and converts lysophosphatidic acid to phosphatidic acid, the second step in de novo phospholipid biosynthesis. Mutations in this gene have been associated with congenital generalized lipodystrophy (CGL), or Berardinelli-Seip syndrome, a disease characterized by a near absence of adipose tissue and severe insulin resistance. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
------	--

基因ID	10555
基因名	AGPAT2
Swiss	O15120
别名	AGPAT2;1-AGPAT2;BSCL;BSCL1;LPAAB;LPAAT-beta

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn