

ANLN Rabbit pAb

货号: B16703

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	120kDa/124kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,MCF7
细胞定位	Cytoplasm,Nucleus,cell cortex,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 845-1124 of human ANLN (NP_061155.2).
序列	ATPLASTSNSLNGDALFTTTFTLQDVSNDFEINIEVYSLVQKKDPSGLDKKKKTSKSKAITPKRLLTSITTKSNIHSSVMASPGGLSAVRTSNFALVGSYTSLSSVGNTKFVLDKVPFLSSLEGHIYLIKICQVNSSVEERGFLTIFEDVSGFGAWHRRWCVLSGNCISYWTPDDEKRKNPIGRINLANCTSRQIEPANREFCARRNTFELITVRPQREDDRETLSQCRDTLCVTKNWLSADTKEERDLWMQKLNQVLVDIRLWQPDACYKPIGKP

靶点信息

研究背景	This gene encodes an actin-binding protein that plays a role in cell growth and migration, and in cytokines is. The encoded protein is thought to regulate actin cytoskeletal dynamics in podocytes, components of the glomerulus. Mutations in this gene are associated with focal segmental glomerulosclerosis 8. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	---

基因ID	54443
基因名	ANLN
Swiss	Q9NQW6
别名	ANLN;FSGS8;Scraps;scra;anillin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn