

MYOZ2 Rabbit pAb

货号: B16717

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	29kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse skeletal muscle,Rat skeletal muscle
细胞定位	Cytoplasm,Z line,myofibril,sarcomere
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-264 of human MYOZ2 (NP_057683.1).
序列	MLSHNTMMKQRKQQATAIMKEVHGNDVDGMDLGKKVSIPRDIMLEELSHLSNRGARLFKMRQRRSDKYTFENFYQSR AQINHSIAMQNGKVDGSNLEGGSSQAPLTPPNTDPDRSPNPDNIAPGYSGPLKEIPPEKFNTTAVPKYYQSPWEQAISND PELLEALYPKLFKPEGKAELPDYRSFNRVATPFGGFEEKASRMVKFKVPDFELLLLTDPRFMSFVNPLSGRRSFNRTPKGWIS ENIPIVITTEPTDDTTVPESEDL

靶点信息

研究背景	The protein encoded by this gene belongs to a family of sarcomeric proteins that bind to calcineurin, a phosphatase involved in calcium-dependent signal transduction in diverse cell types. These family members tether calcineurin to alpha-actinin at the z-line of the sarcomere of cardiac and skeletal muscle cells, and thus they are important for calcineurin signaling. Mutations in this gene cause cardiomyopathy familial hypertrophic type 16, a hereditary heart disorder.
------	---

基因ID	51778
基因名	MYOZ2
Swiss	Q9NPC6
别名	MYOZ2;C4orf5;CMH16;CS-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn