

TFPT Rabbit pAb

货号: B16721

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:100
理论分子量	27kDa/28kDa
实测分子量	27kDa/28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human TFPT (NP_037474.1).
序列	MELEQREGTMAAVGFEEFSAPPGSELALPPLFGGHILESELETEVEFVSGGLGGSGLRERDEEEEAARGRRRRQRELNRRKYQALGRRCREIEQVNERVLNRLHQVQRITRRLQQERRFLMRVLDYSGDDYRASQFTIVLEDEGSQGTDAPTPGNAENEPP EKETLSPPRRTPAPPEPGSPAPGEGPSGRKRRRVPRDGRRAGNALTPELAPVQIKVEEDFGFEADEALDSSWVSRGPKLLPYPTLASPASD

靶点信息

研究背景	Appears to promote apoptosis in a p53/TP53-independent manner.; Putative regulatory component of the chromatin remodeling INO80 complex which is involved in transcriptional regulation, DNA replication and probably DNA repair.
------	---

基因ID	29844
基因名	TFPT
Swiss	P0C1Z6
别名	TFPT;FB1;INO80F;amida

产品验证



Western blot analysis of TFPT expressed in human lung cancer using TFPT Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn