

# C21orf33 Rabbit pAb

货号: B16729

产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB IF/ICC                                                                                              |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000<br><b>IF/ICC:</b> 1:10 - 1:100                                               |
| 理论分子量 | 24kDa/28kDa                                                                                            |
| 实测分子量 | 25KDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | 293T                                                                                                   |
| 细胞定位  | Mitochondrion                                                                                          |
| 纯化    | Affinity purification                                                                                  |

抗原信息


|      |                                                                                                                                                                                                                                      |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 42-268 of human C21orf33 (NP_004640.4).                                                                                                                |
| 序列   | AARVALVLSGCGVYDGTETIHEASAILVHLSRGGAEVQIFAPDVPQMHVIDHTKGQPSEGSRNVLTESARIARGKITDLANLSAANHDAAIFFGGFGAAKNLSTFAVDGKDCKVNKEVERVLKEFHQAGKPIGLCCIAPVLAARKVLRGVEVTVGHEQEEGGKWPYAGTAEAIKALGAKHCVKEVVEAHVDQKNKVTTTAPFMCETALHYIHDGIGAMVRKVLELTGK |

靶点信息

|      |                                                                                                                                                                                                                            |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes a potential mitochondrial protein that is a member of the DJ-1/Pfpl gene family. This protein is overexpressed in fetal Down syndrome brain. Alternate splicing results in multiple transcript variants. |
| 基因ID | 8209                                                                                                                                                                                                                       |

|       |                                         |
|-------|-----------------------------------------|
| 基因名   |                                         |
| Swiss | P30042                                  |
| 别名    | C21orf33;ES1;GATD3;GT335;HES1;KNPH;KNPI |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)