

PIP Rabbit pAb

货号: B16753

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	16kDa
实测分子量	16kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 29-146 of human PIP (N P_002643.1).
序列	QDNTRKIIIKNFDIPKSVRPNDEVTAVLAVQTELKECMVVKTYLISSIPLQGAFNYKYTACLCDDNPKTFYWDFYTNRTVQIA AVVDVIRELGICPDAAVPIKNNRFYTIEILKVE

靶点信息

研究背景	Enables IgG binding activity; aspartic-type endopeptidase activity; and identical protein binding activity. I nvolved in several processes, including detection of chemical stimulus involved in sensory perception of b itter taste; negative regulation of T cell apoptotic process; and proteolysis. Located in extracellular space and nucleus.
基因ID	5304

基因名	PIP
Swiss	P12273
别名	PIP;GCDFP-15;GCDFP15;GPIP4

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn