

C9orf114 Rabbit pAb

货号: B16755

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,K-562
细胞定位	cytoplasm,mitotic spindle,spindle pole centrosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 163-340 of human C9orf114 (NP_057474.2).
序列	HQDLQFAGLLNPLDSPHHMRQDEESEFREGIVVDRPTRPGHGSFVNCGMKKEVKIDKNLEPGLRVTVRLNQQQHPDCKTYHGKVVSSQDPRTKAGLYWGYTVRLASCLSAVFAEAPFQDGYDLTIGTSESGSDVASAQLPNFRHALVVFGGLQGLEAGADADPNLEVAEPSVLFDLYV

靶点信息

研究背景	Required for association of the centrosomes with the poles of the bipolar mitotic spindle during metaphase. Also involved in chromosome alignment. May promote centrosome maturation probably by recruiting A-kinase anchor protein AKAP9 to centrosomes in early mitosis. Binds specifically to miRNA MIR145 hairpin , regulates MIR145 expression at a postranscriptional level.
------	--

基因ID	51490
基因名	SPOUT1
Swiss	Q5T280
别名	C9orf114;CENP-32;HSPC109

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn