

BMX Rabbit pAb

货号: B16773

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	78kDa
实测分子量	78kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	DU145,A375,Jurkat,HeLa,Mouse spleen
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human BMX (N P_001712.1).
序列	MDTKSILEELLKRSQKKKKMSPNNYKERLFVLTKTNLSYHEYDKMKRGSRKGSIEIKKIRCVEKVNLEEQTPVERQYPFQIV YKDGLLYVYASNEESRSQWLKALQKEIRGNPHLLVKY

靶点信息

研究背景	This gene encodes a non-receptor tyrosine kinase belonging to the Tec kinase family. The protein contains a PH-like domain, which mediates membrane targeting by binding to phosphatidylinositol 3,4,5-triphosphate (PIP3), and a SH2 domain that binds to tyrosine-phosphorylated proteins and functions in signal transduction. The protein is implicated in several signal transduction pathways including the Stat pathway, and regulates differentiation and tumorigenicity of several types of cancer cells. Alternatively spliced transcript variants have been found for this gene.
------	---

基因ID	660
基因名	BMX
Swiss	P51813
别名	BMX;ETK;PSCTK2;PSCTK3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn