

BCAP29 Rabbit pAb

货号: B16777

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	28kDa/40kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480,SKOV3,MCF7,THP-1,ES-2,Mouse testis,Mouse heart,Mouse brain,Rat heart
细胞定位	Endoplasmic reticulum membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 125-241 of human BCAP 29 (NP_061332.2).
序列	TQLAKELSNKGVLTQAENTNKAAKKFMEENEKLKRILKSHGKDEECVLEAENKKLVEDQEKLKTELKRTSDALSKAQNDV MEMKMQSERLSKEYDQLLKEHSELQDRLERGNKKRL

靶点信息

研究背景	May play a role in anterograde transport of membrane proteins from the endoplasmic reticulum to the Go lgi. May be involved in CASP8-mediated apoptosis (By similarity.
基因ID	55973

基因名	BCAP29
Swiss	Q9UHQ4
别名	BCAP29;BAP29

产品验证



Western blot analysis of BCAP29 expressed in Mouse testis,Rat heart,SW480 using BCAP29 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn