

SIVA1 Rabbit pAb

货号: B16779

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	11kDa/18kDa
实测分子量	15kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,SW620
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human SIVA1 (NP_068355.1).
序列	MPKRSCPFADVAPLQLKVRVSQRELSRGVCAERYSQEVFDPGVASIACSSCVRAVDGKAVCGQCERALCGQCVRTCW GCGSVACTLCGLVDCSDMYEVLCTSCAMFET

靶点信息

研究背景	This gene encodes a protein with an important role in the apoptotic (programmed cell death) pathway induced by the CD27 antigen, a member of the tumor necrosis factor receptor (TNFR) superfamily. The CD27 antigen cytoplasmic tail binds to the N-terminus of this protein. Two alternatively spliced transcript variants encoding distinct proteins have been described.
基因ID	10572

基因名	SIVA1
Swiss	O15304
别名	SIVA1;CD27BP;SIVA;Siva-1;Siva-2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn