

ZMYND11 Rabbit pAb

货号: B16780

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	60kDa/64kDa/66kDa/70kDa
实测分子量	71kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW620,U-937,Jurkat,Mouse kidney,Mouse lung,Mouse brain
细胞定位	Chromosome,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 278-567 of human ZMYND11 (NP_997644.2).
序列	PNHELWVAKMKGFGFWPAKVMQKEDNQVDVRRFFGHHHQRAWIPSENIQDITVNIHRLHVKRSMGWKKACDELELHQRFLREGRFWKSKNEDRGEEEAESSISSTSNEQLKVTQEPRAKKGRRNQSVEPKKEEPEPETEAVSSSQEIPTMPQPIEKVSVSQTQKKLSASSPRMLHRSTQTTNDGVCQSMCHDKYTKIFNDFKDRMKSDHKRETERVVREALEKLRSEMEEEKRQAVNKA VANMQGEMDRKCKQVKEKCKEEFVEEIKKLATQHKQLISQTKKKQWVNTSLF

靶点信息

研究背景	The protein encoded by this gene was first identified by its ability to bind the adenovirus E1A protein. The protein localizes to the nucleus. It functions as a transcriptional repressor, and expression of E1A inhibits this repression. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	10771

基因名	ZMYND11
Swiss	Q15326
别名	ZMYND11;BRAM1;BS69;MRD30

产品验证



Western blot analysis of ZMYND11 expressed in Mouse kidney,SW620,U-937 using ZMYND11 Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn