

PHF21A Rabbit pAb

货号: B16782

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	70kDa/74kDa
实测分子量	70-80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,293T,Mouse brain,Rat brain
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 461-680 of human PHF2 1A (NP_001095272.1).
序列	TETTTFTFPAPVQPVSLPSPTSTDGDIHEDFCVCRKSGQLLMCDTCSRVIHLDCLDPLKTIPKGMWICPRCQDQMLKKEE AIPWPGTLAIVHSYIAYKAAKEEEKQKLLKWSSDLKQEREQLEQKVQLSNSISKMEMKNTILARQKEMHSSLEKVKQLIR LIHGIDLSKPDSEATVGAISNGPDCTPPANAATSTPAPSPSSQSCTANCNQGEETK

靶点信息

研究背景	The PHF21A gene encodes BHC80, a component of a BRAF35 (MIM 605535)/histone deacetylase (HDAC; s ee MIM 601241) complex (BHC) that mediates repression of neuron-specific genes through the cis-regulat ory element known as repressor element-1 (RE1) or neural restrictive silencer (NRS) (Hakimi et al., 2002 [PubMed 12032298]).[supplied by OMIM, Nov 2010]
基因ID	51317

基因名	PHF21A
Swiss	Q96BD5
别名	PHF21A;BHC80;BM-006

产品验证



Western blot analysis of PHF21A expressed in Mouse brain,Rat brain,U-87MG using PHF21A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn