

TKT Rabbit pAb

货号: B16787

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	67kDa/68kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,293T,Jurkat,MCF7,U-251MG,HepG2,Mouse liver,Mouse kidney,Rat kidney
细胞定位	cytosol,extracellular exosome,nuclear body,nucleoplasm,peroxisome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TKT (NP_001128527.1).
序列	MESYHKPDQQKLQALKDTANRLRISSIQATTAAGSGHPTSCCSAAEIMAVLFFHTMRYKSQDPRNPHNDRFVLSKGHAAP ILYAVWAEAGFLAEAEELNLRKISSDLDGHPVPKQAFSTDVATGSLGQGLGAACGMAYTGKYFDKASYRVYCLLDGELSE GSVWEAMAFASIYKLDNLVAILDINRLGQSDPAPLQHQMIDIYQKRCEAFGWHAIIVDGHSVEELCKAFGQAKHQPTAIIAK TFKGRGITGVEDKESWHGKPLPKNMAEQIIQEISQIQS

靶点信息

研究背景	This gene encodes a thiamine-dependent enzyme which plays a role in the channeling of excess sugar phosphates to glycolysis in the pentose phosphate pathway. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
------	--

基因ID	7086
基因名	TKT
Swiss	P29401
别名	TKT;HEL-S-48;HEL107;SDDHD;TK;TKT1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn