

HSD17B13 Rabbit pAb

货号: B16803

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	29kDa/33kDa
实测分子量	34KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse lung
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 30-211 of human HSD17B13 (NP_835236.2).
序列	RRKSVAGEIVLITGAGHGIGRQTTYEFAKRQSILVLWDINKRGVEETAACRKLGVTAHAYVVDCSNREEIYRSLNQVKKEVG GDTVIVVNNAGTVYPADLLSTKDEEITKTFEVNILGHFWITKALLPSMMERNHGHIVTVASVCGHEGIPYLPYCSSKFAAVG FHRGLTSELQALGKTGI

靶点信息

研究背景	Predicted to enable oxidoreductase activity, acting on the CH-OH group of donors, NAD or NADP as accep tor and steroid dehydrogenase activity. Acts upstream of or within positive regulation of lipid biosynthetic process. Located in lipid droplet.
------	--

基因ID	345275
基因名	HSD17B13
Swiss	Q7Z5P4
别名	HSD17B13;HMFN0376;NIIL497;SCDR9;SDR16C3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn