

AKT1S1/PRAS40 Rabbit pAb

货号: B16814

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100
理论分子量	13kDa/27kDa/29kDa
实测分子量	50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Rat heart
细胞定位	Cytoplasm,cytosol
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-256 of human AKT1S1/PRAS40 (NP_001092103.1).
序列	STDDGSLSEETPAGPPTCSVPPASALPTQQYAKSLPVSPVWGFKEKRTEARSSDEENGPPSSPDLDRIAASMRALVLREA EDTQVFGDLPRPRLNTSDFQKLKRKY

靶点信息

研究背景	AKT1S1 is a proline-rich substrate of AKT (MIM 164730) that binds 14-3-3 protein (see YWHAH, MIM 1135 08) when phosphorylated (Kovacina et al., 2003 [PubMed 12524439]).[supplied by OMIM, Mar 2008]
基因ID	84335

基因名	AKT1S1
Swiss	Q96B36
别名	AKT1S1;Lobe;PRAS40

产品验证



Western blot analysis of AKT1S1/PRAS40 expressed in Mouse heart,Rat heart using AKT1S1/PRAS40 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time : 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn