

TFPI2 Rabbit pAb

货号: B16821

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	25kDa/26kDa
实测分子量	25kDa/26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 23-235 of human TFPI2 (NP_006519.1).
序列	DAAQEPTGNNAEICLLPLDYGPCRALLLRYYDRYTQSCRQFLYGGCEGNANNFYTWACDDACWRIEKVPKVCRLQVSV DDQCEGSTEKYFFNLSSMTCEKFFSGGCHRNRIENRFPDEATCMGFCAPKKIPSCYSPKDEGLCSANVTRYFNPRTCTC DAFTYTGCGGNDNNFVSREDCKRACAKALKKKKKMPKLRFASRIRKIRKKQF

靶点信息

研究背景	This gene encodes a member of the Kunitz-type serine proteinase inhibitor family. The protein can inhibit a variety of serine proteases including factor VIIa/tissue factor, factor Xa, plasmin, trypsin, chymotrypsin and plasma kallikrein. This gene has been identified as a tumor suppressor gene in several types of cancer. Alternative splicing results in multiple transcript variants.
------	--

基因ID	7980
基因名	TFPI2
Swiss	P48307
别名	TFPI2;PP5;REF1;TFPI-2

产品验证



Western blot analysis of TFPI2 expressed in human liver cancer,mouse brain using TFPI2 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blo cking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn