

MLL5 Rabbit pAb

货号: B16840

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	68kDa/98kDa/131kDa/180-204kDa
实测分子量	100KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C2C12
细胞定位	Cell membrane,Cytoplasm,Nucleus speckle,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human MLL5 (NP_891847.1).
序列	MSIVIPLGVDTAETSYLEMAAGSEPESVEASPVVVEKSNSYPHQLYTSSSHSHSHSYIGLPYADHNYGARPPPTPPASPPPSV LISKNEVGIFTTPNFDETSSATTISTSEDGSYGTDVTRCICGFTHDDGYMICCDKCSVWQHIDCMGIDRQHIDPTLYLCERCQ PRNLDKERAVLLQRRKRENMSDGDTS

靶点信息

研究背景	This gene is a member of the myeloid/lymphoid or mixed-lineage leukemia (MLL) family and encodes a protein with an N-terminal PHD zinc finger and a central SET domain. Overexpression of the protein inhibits cell cycle progression. Alternate transcriptional splice variants have been characterized.
------	---

基因ID	55904
基因名	KMT2E
Swiss	Q8IZD2
别名	KMT2E;HDCMC04P;MLL5;NKp44L

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn