

TNFRSF10D Rabbit pAb

货号: **B16842**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	41kDa
实测分子量	47KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 257-386 of human TNFR SF10D (NP_003831.2).
序列	VLFRRRSCPSRVPGAEDNARNETLSNRYLQPTQVSEQEIQQQELAELTGVTVELPEEPQRLLEQAEAEGCQRRRLVPVND ADSADISTLLDASATLEEGHAKETIQDQLVGSEKLFYEDEAGSATSCL

靶点信息

研究背景	The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor contains an extracellular TRAIL-binding domain, a transmembrane domain, and a truncated cytoplamic death domain . This receptor does not induce apoptosis, and has been shown to play an inhibitory role in TRAIL-induced cell apoptosis.
基因ID	8793

基因名	TNFRSF10D
Swiss	Q9UBN6
别名	TNFRSF10D;CD264;DCR2;TRAIL-R4;TRAILR4;TRUNDD

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn