

TDRD3 Rabbit pAb

货号: B16864

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	73kDa/83kDa
实测分子量	73kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U2OS
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human TDRD3 (NP_110421.1).
序列	MLRLQMTDGHISCTAVEFSYMSKISLNTPPGTKVKLSGIVDIKNGFLLNDSNTTVLGGEVEHLIEKWELQRSLSKHNRSNIGTEGGPPPFVFPFGQKCVSHVQVDSRELDRRKTLQVTMPVKPTNDNDEFEKQRTAAIAEVAKSKETKTFGGGGGGARSNLNMNAAGNRNREVLQKEKSTKSEGKHGCVYRELVDEKALKHITEMGFSKEASRQALMDNG

靶点信息

研究背景	Scaffolding protein that specifically recognizes and binds dimethylarginine-containing proteins. In nucleus , acts as a coactivator: recognizes and binds asymmetric dimethylation on the core histone tails associated with transcriptional activation (H3R17me2a and H4R3me2a and recruits proteins at these arginine-methylated loci. In cytoplasm, may play a role in the assembly and/or disassembly of mRNA stress granules and in the regulation of translation of target mRNAs by binding Arg/Gly-rich motifs (GAR in dimethylarginine-containing proteins.
------	---

基因ID	81550
基因名	TDRD3
Swiss	Q9H7E2
别名	TDRD3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn