

NSUN5 Rabbit pAb

货号: B16867

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	36kDa/42kDa/46kDa/50kDa/51kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,RD,Rat brain
细胞定位	nucleolus,nucleoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 127-429 of human NSU N5 (NP_060514.1).
序列	RPGPASQLPRFVRVNTLKTCSDDVVDYFKRQGFSYQGRASSLDDLRLKKGKHFLDPLMPELLVFPAQTDLHEHPLYRAG HLILQDRASCLPAMLLDPPPGSHVIDACAAPGNKTSHLAALLKNQGKIFAFDLDAKRLASMATLLARAGVSCCELAEDFLA VSPSPDRYHEVHYILLDPSCSGSGMPSRQLEEPGAGTSPVRLHALAGFQQRALCHALTFPSLQRLVYSTCSLCQEENEDV VRDALQQNPGAFRLAPALPAWPHRGLSTFPGA EHCLRAS PETTLSSGFFVAVIERVEVPR

靶点信息

研究背景	This gene encodes a member of an evolutionarily conserved family of proteins that may function as methyltransferases. This gene is located in a larger region of chromosome 7 that is deleted in Williams-Beuren syndrome, a multisystem developmental disorder. There are two pseudogenes for this gene located in the same region of chromosome 7. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	--

基因ID	55695
基因名	NSUN5
Swiss	Q96P11
别名	NSUN5;NOL1;NOL1R;NSUN5A;WBSCR20;WBSCR20A;p120;p120(NOL1)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn