

ZCCHC11 Rabbit pAb

货号: B16870

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	81kDa/185kDa
实测分子量	185kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human ZCCHC 11 (NP_001009881.1).
序列	MEESKTLKSENHEPKKNVICEESKAVQVIGNQTLKARNDKSVKEIENSSPNRNSSKKNQNDICIEKTEVKCKVNAANLPGPKDLGLVLRDQSHCKAKKFPNSPVKAEKATISQAKSEKATSLQAKAEKSPKSPNSVKAEKASSYQMKSEKVPSSPAEAEKGPSLLLLKDMRQKTELQQIGKKIPSSFTSVDKVNIEAVGGEKCALQNSPRSQKQQTCTDNTGDSDDDSASGIEDVSDDLKMKNDESNKENSSEMDYLENATVIDESALTPEQRLGLKQAEERLERDHIF

靶点信息

研究背景	ZCCHC11 is an RNA uridylyltransferase (EC 2.7.7.52) that uses UTP to add uridines to the 3-prime end of s ubstrate RNA molecules (Jones et al., 2009 [PubMed 19701194]).
基因ID	23318

基因名	ZCCHC11
Swiss	Q5TAX3
别名	ZCCHC11;PAPD3;TUT4

产品验证



Western blot analysis of ZCCHC11 expressed in 293T using ZCCHC11 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn