

CDKAL1 Rabbit pAb

货号: B16925

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:10 - 1:100
理论分子量	11kDa/54kDa/65kDa
实测分子量	55-65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,HeLa,BT-474,293T,HepG2,Mouse kidney,Mouse heart,Mouse spleen,Mouse lung
细胞定位	Endoplasmic reticulum membrane,Single-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human CDKAL 1 (NP_060244.2).
序列	MPSASCDTLLDDIEDIVSQEDSKPQDRHFVRKDVVPKVRRRNTQKYLQEEENSPPSDSTIPGIQKIWIRTWGCSHNSDG EYMAGQLAAYGYKITENASDADLWLLNSCTVKNPAREDHFRNSIKKAQEENKKIVLAGCVPAQAPRQDYLGKLSIIGVQQID RVVEVVEETIKGHSVRLLGQKKDNGRRLGGARLDLPKIRKNPLIEISINTGCLNACTYCKTKHARGNL

靶点信息

研究背景	The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes.
基因ID	54901

基因名	CDKAL1
Swiss	Q5VV42
别名	CDKAL1

产品验证



Western blot analysis of CDKAL1 expressed in Mouse kidney,HL-60,HeLa using CDKAL1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn