

SMARCB1/SNF5 Rabbit pAb

货号: **B16929**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:500 - 1:1000
理论分子量	43kDa/44kDa
实测分子量	44KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,THP-1,HepG2,HeLa,Mouse testis,Rat testis
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human SMARCB1/SNF5 (NP_003064.2).
序列	LWRRLATVEERKKIVASSHGKTKPNTKDHGYTTLATSVTLKASEVEEILDGNDEKYKAVSISTEPPTYLREQKAKRNSQWVPTLPNSSHHLDVPCSTT

靶点信息

研究背景	The protein encoded by this gene is part of a complex that relieves repressive chromatin structures, allowing the transcriptional machinery to access its targets more effectively. The encoded nuclear protein may also bind to and enhance the DNA joining activity of HIV-1 integrase. This gene has been found to be a tumor suppressor, and mutations in it have been associated with malignant rhabdoid tumors. Alternatively spliced transcript variants have been found for this gene.
------	--

基因ID	6598
基因名	SMARCB1
Swiss	Q12824
别名	SMARCB1;BAF47;CSS3;INI1;MRD15;PPP1R144;RDT;RTPS1;SNF5;SNF5L1;SWNTS1;Sfh1p;Snr1;hSNFS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn