

FUT2 Rabbit pAb

货号: B16950

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:1000
理论分子量	37kDa/39kDa
实测分子量	35KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29
细胞定位	Golgi apparatus,Golgi stack membrane,Single-pass type II membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 29-150 of human FUT2 (NP_001091107.1).
序列	QQRLAKIQAMWELPVQIPVLASTSKALGPSQLRGMWTTINAIGRLGNQMGEYATLYALAKMNGRPAFIPAQMHSSTLAPIFRITLPVLHSATASRIPWQNYHLNDWMEEEYRHIPGEYVRFTGY

靶点信息

研究背景	The protein encoded by this gene is a Golgi stack membrane protein that is involved in the creation of a p precursor of the H antigen, which is required for the final step in the soluble A and B antigen synthesis path way. This gene is one of two encoding the galactoside 2-L-fucosyltransferase enzyme. Two transcript vari ants encoding the same protein have been found for this gene.
基因ID	2524

基因名	FUT2
Swiss	Q10981
别名	FUT2;B12QTL1;SE;SEC2;Se2;sej

产品验证



Western blot analysis of FUT2 expressed in HT-29 using FUT2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn