

TP53BP2 Rabbit pAb

货号: B16958

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	IHC: 1:50 - 1:200
理论分子量	111kDa/125kDa/126kDa
实测分子量	111kDa/125kDa/126kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,Nucleus,perinuclear region
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 360-600 of human TP53 BP2 (NP_001026855.2).
序列	STMPRMPSRPELLVKPALPDGSLVIQASEGPMKIQTLPNMRSGAASQTKGSKIHPVGPDWSPSNADLFPSQGSASVPQST GNALDQVDDGEVPLREKEKKVRPFMFDAVDQSNAPPSFGTLRKNQSSDILRDAQVANKNVAKVPPPVP TKPKQINLPY FGQTNQPPSDIKPDGSSQQLSTVVPSMGTKPKPAGQQPRVLLSPSIPSVGQDQTLSPGSKQESPPAAAVRPFTPQPSKDT L

靶点信息

研究背景	This gene encodes a member of the ASPP (apoptosis-stimulating protein of p53) family of p53 interacting proteins. The protein contains four ankyrin repeats and an SH3 domain involved in protein-protein interactions. It is localized to the perinuclear region of the cytoplasm, and regulates apoptosis and cell growth through interactions with other regulatory molecules including members of the p53 family. Multiple transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	7159
基因名	TP53BP2
Swiss	Q13625
别名	TP53BP2;53BP2;ASPP2;BBP;P53BP2;PPP1R13A

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn