

IQGAP1 Rabbit pAb

货号: B16966

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:100
理论分子量	189kDa
实测分子量	189kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,293T,HepG2,HT-29,HeLa
细胞定位	Cell membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 500-750 of human IQGA P1 (NP_003861.1).
序列	ELMKLKAQAHAENNEFITWNDIQACVDHVNLVVQEEHERILAIGLINEALDEGDAQKTLQALQIPAAKLEGVLAEVAQHYQD TLIRAKREKAQEIQDESAVLWLDEIQGGIWQSNKDTQEAQKFALGIFAINAEVESGDVGKTLSALRSPDVGLYGVIEPCGET YHSDLAEAKKKKLAVGDNNSKWVKHWVKGGYYYYHNLETQEGGWDEPPNFVQNSMQLSREEIQSSISGVTAAYNREQL WLANEGLIT

靶点信息

研究背景	This gene encodes a member of the IQGAP family. The protein contains four IQ domains, one calponin ho mology domain, one Ras-GAP domain and one WW domain. It interacts with components of the cytoskele ton, with cell adhesion molecules, and with several signaling molecules to regulate cell morphology and motility. Expression of the protein is upregulated by gene amplification in two gastric cancer cell lines.
------	---

基因ID	8826
基因名	IQGAP1
Swiss	P46940
别名	IQGAP1;HUMORFA01;SAR1;p195

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn