

# IL3 Rabbit pAb

货号: B16969

## 产品信息

|       |                                                                                                      |
|-------|------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse                                                                                          |
| 宿主    | Rabbit                                                                                               |
| 克隆性   | Polyclonal                                                                                           |
| 预测反应  |                                                                                                      |
| 应用    | WB                                                                                                   |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:1000                                                                            |
| 理论分子量 | 17kDa                                                                                                |
| 实测分子量 | 17kDa                                                                                                |
| 形式    | Liquid                                                                                               |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                         |
| 阳性对照  | Raji,Mouse thymus                                                                                    |
| 细胞定位  | Secreted                                                                                             |
| 纯化    | Affinity purification                                                                                |

## 抗原信息


|      |                                                                                                                                    |
|------|------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 20-152 of human IL3 (NP_000579.2).                   |
| 序列   | APMTQTTPPKTSWVNCSNMIDEIITHLKQPPLPLDFNNLNGEDQDILMENNLRPNLEAFNRAVKSLQNASAIESILKNLLPCLPLATAAPTRHPIHIKDGDWNEFRRLTFYKLTLENAQAQQTTLSLAIF |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                                           |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is a potent growth promoting cytokine. This cytokine is capable of supporting the proliferation of a broad range of hematopoietic cell types. It is involved in a variety of cell activities such as cell growth, differentiation and apoptosis. This cytokine has been shown to also possess neurotrophic activity, and it may be associated with neurologic disorders. |
| 基因ID | 3562                                                                                                                                                                                                                                                                                                                                                                                                      |

|       |                         |
|-------|-------------------------|
| 基因名   | IL3                     |
| Swiss | P08700                  |
| 别名    | IL3;IL-3;MCGF;MULTI-CSF |

产品验证



Western blot analysis of IL3 expressed in Raji,Mouse thymus using IL3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)