

OCA2 Rabbit pAb

货号: B16972

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	74kDa/90kDa/92kDa
实测分子量	65-80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-937,Mouse eye
细胞定位	Melanosome membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human OCA2 (NP_000266.2).
序列	MHLEGRDGRRYPGAPAVELLQTSVPSGLAELVAGKRRLPRGAGGADPSHSCPRGAAGQSSWAPAGQEFASFLTGRSH SSLPQMSSSRSKDSCFTENTPLLRNSLQEKGSRCIPVYHPEFITAEESWEDSSADWERRYLLSREVSGLSASASSEKGDLL DSPHIRLRLSKLRRVCQWLKV

靶点信息

研究背景	This gene encodes the human homolog of the mouse p (pink-eyed dilution) gene. The encoded protein is believed to be an integral membrane protein involved in small molecule transport, specifically tyrosine, which is a precursor to melanin synthesis. It is involved in mammalian pigmentation, where it may control skin color variation and act as a determinant of brown or blue eye color. Mutations in this gene result in type 2 oculocutaneous albinism. Alternative splicing results in multiple transcript variants.
------	--

基因ID	4948
基因名	OCA2
Swiss	Q04671
别名	OCA2;BEY;BEY1;BEY2;BOCA;D15S12;EYCL;EYCL2;EYCL3;HCL3;P;PED;SHEP1;P protein

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn