

SH2B1 Rabbit pAb

货号: B16996

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC IF/ICC
推荐浓度	IHC: 1:50 - 1:200 IF/ICC: 1:20 - 1:50
理论分子量	70kDa/72kDa/79kDa
实测分子量	70kDa/72kDa/79kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm,Membrane,Nucleus
纯化	Affinity purification

抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 360-520 of human SH2B 1 (NP_056318.2).
序列	DAQHVKAWSVDIQECLSPGPCPATSPRPMTLPLAPGTSFLTRENTDSLELSCLNHSESLPSQDLLLGPSESNDRLSQGAY GGLSDRPSASISPSSASIAASHFDSMELLPPELPPRIPIEEGPPTGTVHPLSAPYPPLDTPETATGSFLFQGEPEGGEDQ

靶点信息

研究背景	This gene encodes a member of the SH2-domain containing mediators family. The encoded protein media tes activation of various kinases and may function in cytokine and growth factor receptor signaling and ce llular transformation. Alternatively spliced transcript variants have been described.
基因ID	25970

基因名	SH2B1
Swiss	Q9NRF2
别名	SH2B1;PSM;SH2B

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn