

TFCP2 Rabbit pAb

货号: B17001

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	51kDa/57kDa
实测分子量	64kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,THP-1,SKOV3,K-562,Mouse kidney,Mouse brain,Mouse testis
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TFCP2 (NP_005644.2).
序列	MAWALKLPLADEVIESGLVQDFDASLSGIGQELGAGAYSMSDVLALPIFKQEESLPPDNENKILPFQYVLCAATSPAVKLH DETLYLNQGQSYEIRMLDNRKLGELPEINGKLVKSIFRVVFHDRRLQYTEHQQLLEGWRWNRPGDRILDIDIPMSVGIIDPR ANPTQLNTVEFLWDPAKRTSVFIQVHCISTEFTMRKHGGEKGVFPFRVQIDTFKENENGEYTEHLHSASCQIKVFKPKGADR KQKTDREKMEKRTPHEKEYQPSYETTILTECSPW

靶点信息

研究背景	This gene encodes a transcription factor that binds the alpha-globin promoter and activates transcription of the alpha-globin gene. The encoded protein regulates erythroid gene expression, plays a role in the transcriptional switch of globin gene promoters, and it activates many other cellular and viral gene promoters. The gene product interacts with certain inflammatory response factors, and polymorphisms of this gene may be involved in the pathogenesis of Alzheimer's disease.
------	--

基因ID	7024
基因名	TFCP2
Swiss	Q12800
别名	TFCP2;LBP1C;LSF;LSF1D;SEF;TFCP2C

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn