

SCNN1A Rabbit pAb

货号: B17009

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	28kDa/73kDa/75kDa/77kDa/78kDa/81kDa
实测分子量	75KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,A-549
细胞定位	Apical cell membrane,Cell projection,Multi-pass membrane protein,cilium
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 105-225 of human SCN N1A (NP_001029.1).
序列	LFGEYFSYPVSLNINLNSDKLVFPAVTICTLNPYRYPEIKEELEELDRITEQTLFDLYKYSSFTTLVAGSRSRRLRGTLPHPLQLRLRVPPPPHGARRARSVASSLRDNNPQVDWKDWKIG

靶点信息

研究背景	Nonvoltage-gated, amiloride-sensitive, sodium channels control fluid and electrolyte transport across epit helia in many organs. These channels are heteromeric complexes consisting of 3 subunits: alpha, beta, a nd gamma. This gene encodes the alpha subunit, and mutations in this gene have been associated with p pseudohypoaldosteronism type 1 (PHA1), a rare salt wasting disease resulting from target organ unrespon siveness to mineralocorticoids. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
------	---

基因ID	6337
基因名	SCNN1A
Swiss	P37088
别名	SCNN1A;BESC2;ENaCa;ENaCalpha;SCNEA;SCNN1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn