

UMPS Rabbit pAb

货号: B17019

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	22kDa/33kDa/42kDa/52kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat
细胞定位	cytoplasm,cytosol,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 181-480 of human UMP S (NP_000364.1).
序列	EILEQQKKVDAETVGRVKRFIQENVFVAANHNGSPLSIKEAPKELSFGARAEIPRIHPVASKLLRLMQKKETNLCLSADVSLARELLQLADALGPSICMLKTHVDILNDFTLDVMKELITLAKCHEFLIFEDRKFDIGNTVKKQYEGGIFKIASWADLVNAHVVP GSGVVKGLQEVLPLHRGCLLIAEMSSTGSLATGDYTRAAVRMAEEHSEFVVGFIGSRVSMKPEFLHILTPGVQLEAGGD NLGQQYNPQEVIGKRGSDIIIVGRGIISAADRLEAAEMYRKAWEAYLSRLGV

靶点信息

研究背景	This gene encodes a uridine 5'-monophosphate synthase. The encoded protein is a bifunctional enzyme that catalyzes the final two steps of the de novo pyrimidine biosynthetic pathway. The first reaction is carried out by the N-terminal enzyme orotate phosphoribosyltransferase which converts orotic acid to orotidine-5'-monophosphate. The terminal reaction is carried out by the C-terminal enzyme OMP decarboxylase which converts orotidine-5'-monophosphate to uridine monophosphate. Defects in this gene are the cause of hereditary orotic aciduria. Alternate splicing results in multiple transcript variants.
------	---

基因ID	7372
基因名	UMPS
Swiss	P11172
别名	UMPS;OPRT

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn