

PP2A-B56γ/PR61γ/PPP2R5C Rabbit pAb

货号: B17021

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	52kDa/56kDa/61kDa/62kDa/64kDa
实测分子量	68-71kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,K-562,Mouse testis
细胞定位	Chromosome,Nucleus,centromere
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human PP2A-B56γ/PR61γ/PP2A-B56γ/PR61γ/PPP2R5C (NP_848702.1).
序列	MLTCNKAGSRMVVDAANSNGPFQPVVLLHIRDVPPADQEKLFIQKLRQCCVLFDFVSDPLSDLKWKEVKRAALSEMVEYI THNRNVITEPIYPEVVHMFVNMFRSLPPSSNPTGAEFDPEEDEPTLEAAWPHLQLVYEFFFLRFLESPPDFQPNIKKYIDQKF VLQELLEFDSEDPREDFLKTTLHRIYGKFLGLRAYIRKQINNIFYRFIYETEHNGIAELLEILGSIINGFALPLKEEHKIFLLKVL LPLHKVKSLSVYHPQLAYCVVQFLEKDSTLTPVVMALLKYWPKTHSPKEVMFLNELEEILDVIEPSEF

靶点信息

研究背景	The product of this gene belongs to the phosphatase 2A regulatory subunit B family. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. This gene encodes a gamma isoform of the regulatory subunit B56 subfamily. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	5527
基因名	PPP2R5C
Swiss	Q13362
别名	PPP2R5C;B56G;PR61G;B56gamma

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn