

GGH Rabbit pAb

货号: B17028

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	35kDa
实测分子量	32kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse kidney
细胞定位	Lysosome,Melanosome,Secreted,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 25-230 of human GGH (NP_003869.1).
序列	RPHGDTAKKPIIGILMQKCRNKVMKNYGRYYIAASYVKYLESAGARVVPVRLDLTEKDYEILFKSINGILFPGGSVDLRRSDY AKVAKIFYNLSIQSFDDGDYFPVWGTCLGFEELSLISGECLLTATD TVDVAMPLNFTGGQLHSRMFQNFPTELL LSLAVEP LTANFHKWSLSVKNFTMNEK LKKFFNVLTTNTDGKIEFIST

靶点信息

研究背景	This gene catalyzes the hydrolysis of folylpoly-gamma-glutamates and antifolylpoly-gamma-glutamates b y the removal of gamma-linked polyglutamates and glutamate.
基因ID	8836

基因名	GGH
Swiss	Q92820
别名	GGH;GH

产品验证



Western blot analysis of GGH expressed in HeLa,Mouse kidney using GGH Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn